HLS MODULE PHYSIOLOGY LAB 1 HEMATOLOGICAL TESTS

BY

Dr. Fatma Farrag Ali
Associate Professor of Medical Physiology
Faculty of Medicine-Mutah University
2024-2025



1- Determination of Hemoglobin (Hb) Content

Apparatus: Sahli hemoglobinometer or haemometer.


Principle:

This methods depends on measurement of the Hb content (or its % of normal) by matching the blood color with a standard solution that has the color of normal blood (containing 15 gm Hb%).

Sahli apparatus consists of:

- 1. A pipette to measure 20 cu mm (mm³) OR (20 μl) of blood.
- 2. A central graduated diluting tube; two graduations are present:
- Yellow colored graduation indicating amount of Hb in g/100 ml blood.
- Red colored graduation indicating Hb% of normal in blood (the normal is 100% = 15 gm/100 ml blood).
- 3. Two lateral colored standard tubes.







Method:

- 1. 0.1N (N/10) HCl is added to the diluting tube up to the mark 10%.
- 2. Sterilize the finger with alcohol; allow drying and making a pin prick.
- 3. Suck blood with the pipette up to the mark **20 cu mm (0.02 ml; 20** µl).
- 4. Evacuates the pipette into the diluting tube to which 0.1N HCL was previously added. **HCl converts Hb to acid hematin (brownish)**.
- 5. Distilled water is then added drop by drop to the sample till its color becomes the same as that of the standard tubes.
- 6. The Hb content (in gm) and its % of normal can then be directly obtained from gradations on the dilution tube.

RESULTS

The amount of hemoglobin in:

- 1. Normal adult male: is 14-17 gm %. (average: 15 gm %)
- 2. Normal adult female: is 12-15.5 gm%. (average: 13.5 gm %)
- 3. Newborn: is 19-20 gm%.

Clinical significance:

Diagnosis of anemia or polycythemia.

2- Measurement of Hematocrit (Hct) Or Packed Cell Volume (PCV)



Hematocrit (Hct)


Definition: Percentage ratio of RBCs volume to total blood volume.


OR the volume of packed cells; RBCs (PCV) in 100 ml Blood.

So, it is also called **PCV**.

Hct (PCV) =
$$\frac{\text{Packed cell volume}}{\text{Total blood volume}} \times 100$$


Method of estimation: graduated heparinized capillary tube.





Principle:

- Separation of blood cells from plasma by centrifugation depending on differences in specific gravity.
- RBCs have specific gravity 1090 will sink into the bottom with WBCs and platelets (specific gravity 1060) and plasma (1025-1030) on top.




Steps of measurement of Hct (PCV)

- A blood sample from the subject is placed in a graduated heparinized capillary tube and is then centrifuged.
- The red cells are packed in the bottom of the tube while the clear plasma remains above and the white cells and platelets form a small buffy layer just above the red cell column.
- Hct is then calculated by dividing the red blood cells column by the total blood column and multiplying by 100.



Hematocrit Determination



Where:

- H₁ = height of the RBC column
- H₂ = height of the RBC
 + height of the plasma column
- Calculate Hc% (hematocrit) value

Calculation:

RESULTS

Normally, Hct is about :

- ■45% for adult male.
- ■35% for adult female.
- ■60% for newborn (due to polycythemia).

Changes in Hct


Hct is affected by changes in red cell volume relative to the plasma volume.

Decreases in:

- 1- Anemia
- 2- Hydration
- 3- After fluid transfusion

Increases in:

- 1- Polycythemia
- 2- Dehydration
- 3- Loss of plasma as in burns



3- Measurement of Erythrocyte Sedimentation Rate (ESR)




Definition

It is the distance sedimented by R.B.Cs in **mm** when put in a vertical stationary tube at the end of 1 hour.


Method of ESR determination


Westergren tube



WESTERGREN TUBE


- ✓ The Westergren tube is open at both ends.
- ✓ It is 30 cm in length and 2.5 mm in diameter.
- ✓ The lower 20 cm are marked with 0 at the top and 20 at the bottom.






Principle:

- ✓ When anticoagulated blood is allowed to stand in a narrow vertical glass tube, undisturbed for a period of time, the RBCs settle out from the plasma.
- ✓ The rate at which they settle is measured as the number of millimeters (mm) of clear plasma present at the top of the column after one hour (mm/hr).
- ✓ RBCs tend to sediment towards the bottom because of:
- 1. RBCs density is greater than that of plasma.
- 2. RBCs tend to aggregate forming rouleaux shapes in which the ratio of mass to surface area increases and being heavier than plasma, they sediment to the bottom more rapidly.



REQUIREMENTS:

- Anticoagulated blood (0.5 ml of 3.8% sodium citrate solution + 2 ml blood): 1 part citrate to 4 parts of blood (anticoagulant / blood ratio: 1:4).
- Westergren tube.
- Westergren stand.
- Rubber bulb (sucker).




Procedure:

- Mix the anticoagulated blood.
- Draw the blood into Westergren tube up to 0 mark with the help of rubber bulb.
- Set the tube upright (vertical position) in the stand.
- Leave the tube undisturbed, the RBCs settle down leaving a clear plasma column above.
- At the end of 1 hour, read the result.



Normal values of ESR:

Sex 1st hour (average values)
Male 5 mm/hr
Female 8 mm/hr



Factors that affect ESR:

Factors that increase ESR:

- 1. Physiological: as during pregnancy and menstruation, and after meals.
- **2. Pathological:** ESR generally increases in cases of **inflammation and tissue destruction** due to increased fibrinogen or globulins e.g.
- Acute & chronic infection: As tonsillitis, appendicitis and tuberculosis.
- Myocardial infarction.
- Degenerative & neoplastic diseases
- Rheumatoid arthritis.

Factors that decrease ESR:

Pathological decrease in: polycythemia.

N.B.: In anemia, ESR increases except in iron deficiency anemia and spherocytosis (ESR decreases due to reduction of intrinsic ability of RBCs to sediment).

Clinical importance of ESR estimation:

The erythrocyte sedimentation rate (ESR) is a non-specific test. WHY?

It is raised in a wide range of infectious, inflammatory, degenerative, and malignant conditions associated with changes in plasma proteins, particularly increases in fibrinogen and immunoglobulins.

ESR is not a diagnostic test of any particular disease but used to indicate:

- ✓ The presence of tissue damage.
- ✓ The severity of the disease (ESR is increased above 100 mm/hr in rheumatic fever, tuberculosis (T.B.) and malignancy where tissue damage is high).
- ✓ The follow up (the activity of disease and its response to treatment (prognostic test).

