



URINE ANALYSIS
URINARY TRACT INFECTIONS
UGT MODULE LAB 1
2023-2024

Mathhar Ahmad Abu Morad M.A.Sc.

Routine Urine Culture

Urine:

- Urine carries waste products and excess water out of the body.
- Normal urine is typically pale yellow and clear.
- Obvious abnormalities in the color, clarity, and cloudiness may suggest different diseases.



Normal Urine



Abnormal Urine

Routine Urine Culture

Aim of the test

- An etiological diagnosis of bacterial urinary tract infection with identification and susceptibility test of the isolated bacteria(s).


Types of specimen

- Urine (Midstream urine), suprapubic aspiration, catheterized urine.
- ✗ **Note:** First morning specimens yield highest bacterial counts from overnight incubation in the bladder and are the best specimens.

Criteria of specimen rejection

- Un-refrigerated specimen older than **2 hours** may be subject to overgrowth and may not yield valid results; unlabeled specimen; mislabeled specimen; specimen in expired transport container; 24 hours urine specimens.

Urine Analysis



Specimen Collection

Patient

-Collecting **preparing** of midstream urine for investigation:

☐ Patient not needing assistance:

✘ Give the patient a suitable container.

✘ Instruct the patient to collect the midstream urine .

✘ Tell the patient not to touch the inside or rim of the container.

✘ Tell the patient to close the container properly.




Specimen Collection

Who will collect the specimen

- Midstream urine is collected by the patient.
- If disabled, nursing staff will assist in collection.
- For catheterized specimen, nursing staff will collect the specimen.
- Suprapubic aspiration is performed by the physician.

Quantity of specimen

To fill line in transport tube (~20 mL).



Suprapubic aspiration

Transportation

Time relapse before processing the sample

The maximum time allowed for processing a urine sample is **2 hours** from the time of collection.

Storage

At room temperature unless delay is inevitable; it must be refrigerated or mixed with preservative like boric acid.

Processing- Nonculture Methods

Macroscopic Urinalysis

Macroscopic examination used to view elements that are visible by naked eye.

1- Hematuria: is the presence of abnormal numbers of red cells in urine due to:

- a. Glomerular damage.
- b. Tumors.
- c. Urinary tract stones.
- d. Upper and lower urinary tract infections.

Processing- Nonculture Methods

Macroscopic Urinalysis

Hematuria

Two Types of Hematuria

- **Gross hematuria:** means that the blood can be seen by the naked eye. The urine may look pinkish, brownish, or bright red.



**Gross
Hematuria**

Processing- Nonculture Methods

Macroscopic Urinalysis

2- Hemoglobinuria:

- Presence of hemoglobin in urine due to rupturing of RBCs
- This may occur in malaria, typhoid, yellow fever, hemolytic jaundice and other diseases.



Processing- Nonculture Methods

Macroscopic Urinalysis

3- Pyuria :

Refers to the presence of abnormal numbers of leukocytes that may appear with infection in either the upper or lower urinary tract or with acute glomerulonephritis.



Processing- Nonculture Methods

Microscopic Urinalysis

The primary purpose of microscopic examination of urine sediment is to detect abnormal formed elements (eg, cells, casts, crystals) in the sample.

Microscopic hematuria & pyuria

Microscopic hematuria & pyuria means that the urine is clear, but RBCs and WBCs can be seen only under a microscope.

Pyuria: refers to urine which contains pus cells granulocytes.

Normal values:

- Men: <2 WBCs per high power field
- Women: <5

Normal values for RBCs in urine:

4 RBCs per high power field (RBC/HPF).



Microscopic Hematuria

Processing- Nonculture Methods

Microscopic Urinalysis

Microscopic hematuria & pyuria

Reporting:

WBC count recorded as:

<10/ml

10-100/ml

100-500/ml

>500/ml

Other findings may be recorded as:




-/+ = Scanty

+ = Few

++ = Moderate

+++ = Many

WBCs and RBCs in Urine



RBCs in urine

These white blood cells in urine have lobed nuclei and refractile cytoplasmic granules.

Processing- Nonculture Methods

Microscopic Urinalysis

Bacteria

- Bacteria are common in urine specimens (from contamination).
- Therefore, micorganisms in carefully collected urines should be interpreted in view of clinical symptoms.



Microscopic Urinalysis


Urine Crystals

Importance of the urine crystals (crystalluria):

1. These crystals are important in the case of kidney stones.
2. Renal damage was caused by the crystals.
3. In liver diseases.
4. Inborn error of metabolism.

Reporting of the crystalluria:

1. Rare/HPF.
2. FeW/HPF.
3. Moderate/HPF.
4. Many/HPF.



Microscopic Urinalysis

Casts




Are formed only in the **distal convoluted tubule (DCT)** or the **collecting duct**

The major component is the Tamm-Horsfall protein
Other proteins are albumin and immunoglobulins.


Cellular casts are made from protein and contain white blood cells, red blood cells, or epithelial cells. **Non-cellular casts** are made from protein and can contain fat. Hyaline, fatty, granular, and waxy casts are all non-cellular casts.

Cellular casts, fatty casts, granular casts, and waxy casts are not normally present in urine, so their presence could indicate the patient has kidney problems.

Hyaline casts are made from only protein and can typically be found in low numbers in urine, with 0 to 4 casts per high-powered field. Larger numbers of hyaline casts could indicate dehydration, physical exertion, fever, or kidney disease



UTI testing pathway



Processing- Nonculture Methods

Detection of bacteriuria by nitrite test

- Used for screening for **bacteria**.
- Normal urine contain nitrate but not nitrites.
- In the presence of bacteria, the normally present nitrate in the urine is reduced to nitrite.




- Positive test indicates presence of more than 10 organisms/ml.
- Detected by **dipstick chemical analysis**

Processing- Nonculture Methods

Detection of pyuria by leukocyte esterase tests

-Depends on esterase method:



+ve result: means more than 5 leucocytes/hpf. (high power field)

-Detected by dipstick chemical analysis

Dipstick chemical analysis

The squares on the dipstick represent the following components in the urine

	Glucose
	Bilirubin
	Ketones
	Specific Gravity
	Blood
	pH
	Protein
	Urobilinogen
	Nitrite
	Leukocyte Esterase



Dipstick chemical analysis

- Glucose
- Bilirubin
- Ketone
- Specific Gravity
- Blood
- Protein
- Urobilinogen
- Nitrite
- Leukocyte
- pH

Glucose 30 seconds	Negative	g/dl (%)	1/10 (tr.)	1/4	1/2	1	>=2
		mg/dl	100	250	500	1000	>=2000
Bilirubin 30 seconds	Negative		Small +	Moderate ++	Large +++		
Ketone 40 seconds	Negative	mg/dl	trace 5	small 15	mod. 40	large 80	large 160
Sp. Gr. 45 seconds	1000	1005	1010	1015	1020	1025	1030
Blood 60 seconds	Negative	Non-hemo-lyzed trace	Non-hemo-lyzed mod.	hemo-lyzed trace	small +	mod. ++	Large +++
Protein 60 seconds	Negative	trace	mg/dl 30	100	300	>=2000	
			+	++	+++	++++	
Urobilinogen 60 seconds	Normal 0.2	Normal 1	mg/dl 2	4	8	(1mg = approx. 1EU)	
Nitrite 60 seconds	Negative	Positive	Positive	(Any degree of uniform pink colour is positive)			
Leukocytes 2 minutes	Negative	trace	small +	mod. ++	Large +++		
pH 60 seconds	5.0	6.0	6.5	7.0	7.5	8.0	8.5

Dipstick chemical analysis

leukocyte esterase tests




Leukocytes: Indicates infection or inflammation

Normal=negative

- **Pyuria:** Leukocytes in urine
- **Cystitis:** Bladder infection
- **Pyelonephritis:** Kidney infection

Dipstick chemical analysis

Nitrite test




Normal=negative

Dipstick chemical analysis

Dipstick Urinalysis Interpretation- Blood

Blood: Almost always indicates pathology because RBC are too large to pass through glomerulus



Normal=negative

- Hematuria: Blood in urine
- Possible causes: Kidney stone, infection, tumor
- **Caution:** Very common finding in women because of menstruation.

الخدمات الطبية الملكية

مستشفى : _____

فحص البول URINALYSIS

اسم المريض الكامل : _____

الرقم الطبي : _____ الخلية : _____

الرقم الوطني :

الرتبة : _____


العمر : _____

الجنس : ذكر انثى

الزواج : متزوج أعزب


الطبيب الامتصاصي :	اسم و توقيع الطبيب المشرف :	القسم / العيادة :	تاريخ الطلب :
Diagnosis & Relevant Information :			
Test	Result	Test	Result
Albumin	+ (1)	Microscopic	
Sugar	+ + (2)	sediment :	
Other Tests :		R. B. C	8-10
		W. B. C	8-20
		casts	
		Crystals: <i>Ca²⁺</i>	<i>oxalate</i>
		Other Findings	
Comments :			
اسم و توقيع الطبيب المختبر :	اسم و توقيع المشرف :	التاريخ :	التسلسل :

Urine Analysis






Urine Analysis

Culture Methods



Culture of urine sample



General Criteria to Diagnose UTI

Suprapubic Aspiration:

Any growth.

Catheterization:

$\geq 10,000$ colony forming units/ml.

Midstream Clean Catch:

$\geq 100,000$ colony forming units/ml.

Urine Analysis

Culture Methods

- Significant Growth of E. coli in MacConkey Agar



الخدمات الطبية الملكية

مستشفى :

فحص الأحياء الدقيقة
MICROBIOLOGY

اسم المريض الكامل :

الرقم الطبي : الفئة :

الرقم الوطني :

الرتبة : أنثى ذكر

العمر : أعزب متزوج

تاريخ الطلب:	القسم/العيادة:	اسم وتوقيع الطبيب المشرف:	الطبيب الاختصاصي:
Diagnosis :		Type & Source of Specimen:	Test Requested:
Sensitivity Test		For Lab. Use	
Sensitive	Resistant	Result : E. coli > 10 ⁵	
Aug	Am		
Ta ₂	GN		
AK	NA		
CP	SXT		
oFX			
اسم وتوقيع طبيب المختبر:	اسم وتوقيع المرفق:	التاريخ:	التسلسل:

Pathogens and commensals

Common pathogens	Commensal flora
<i>Neisseria gonorrhoeae</i> any colony on chocolate or TM agar (special request).	<i>Diphtheroid bacilli</i>
<i>E.coli</i> and other <i>Enterobacteriaceae</i>	<i>Lactobacillus spp</i>
<i>Enterococcus spp</i>	<i>Coagulase negative Staphylococci</i>
* <i>Staphylococcus aureus</i> Pure culture regardless to the no. of CFUs.	<i>Alpha Haemolytic Streptococci</i>
<i>Staph saprophyticus</i>	<i>Bacillus spp</i>
<i>Corynebacterium jeikeium</i>	<i>Non pathogenic Neisseria spp.</i>
<i>Acinetobacter spp</i>	<i>Anaerobic cocci</i>
<i>Pseudomonas spp</i>	<i>Commensal Mycobacterium</i>
* <i>Gardnerella vaginalis</i> Unusual	<i>Commensal Mycoplasma spp.</i>
<i>Beta -haemolytic streptococci</i>	
* <i>Salmonella spp</i> (early stage of infection)	* yeast
Parasites	
<i>Schistosoma haematobium</i>	
<i>Trichomonas vaginalis</i>	

*** Diagnostic Microbiology, BAILEY & SCOTT, 9th
EDITION**